
ComPact ecoPIN

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ComPact eco PIN – your tool for fast, non-destructive and accurate Karat-measurement of Precious Metals and finished Jewellery.

Based on over 20-years of experience in non-destructive material analysis the *ComPact eco PIN* was developed. Usage of a high-resolution silicon detector in conjunction with digital signal processing the *ComPact eco PIN* achieves best detection limits and highest precision. Based on an about 4 times better resolution (in comparison with conventional technique) the element separation is much better. Furthermore lowest detection limits are achieved because of better signal to noise ratio.

The *ComPact eco PIN*, equipped with modern, user friendly software and exact evaluation algorithm, is the proper analyser for precisely calculating the composition of solid pieces of jewellery, precious coins and other noble metals. Simply position the piece of interest with the aid of a video microscope and start the measurement. The analysis is performed non-contact and non-destructive in about 2 to 3 minutes. The result is given in weight-% or Karat and might be printed as a certificate.

The usage is easy and safe, even for untrained sales personnel at jewellery showrooms. A special laboratory room is not needed.

ComPact eco PIN - fast return on investment!

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What accuracy can YOU expect on gold alloys?

Calibration ranges

The *ComPact eco PIN* is calibrated for gold, silver and platinum. For calibration certified reference samples are used. Calibration is valid for the concentration ranges shown in table no. 1.

Alloy system	Au	Ag	Cu	Zn	Pd	Ni
Gold	30 - 100	0 - 70	0 - 40	0 - 10	0 - 30	0 - 20

Table 1: Cal ibrated co ncentration ranges

Accuracy of analysis

An accuracy for Au of $\pm 0.1\%$ -or 0.05 Karat- can be achieved with an analysis time of 2 minutes. For all other admixtures the accuracy is in the same range. The accuracy of the X-ray spectrometric determination of gold shall be validated by certified standards of jewellery alloy samles. The following table 2 enables to compare X-ray measured values (is) with these of gold standards (given).

Alloy system	Au		Ag		Cu	
	given	is	given	is	given	is
Au1	33,32	32,40	66,59	67,59		
Au2	33,35	33,38	44,65	45,19	21,98	21,42
Au13	91,67	91,67	2,76	2,98	5,28	5,35
Au15	98,60	98,78			1,40	1,21
Au16	99,99	100,0				

Table 2: Analysis results

Technical Data

Tube power	max. 50 Watt, application optimised	HV-Generator:	max. 50kV, 1.2mA, programmable
X-ray tube:	long life, high efficiency tungsten target, air cooled optional motorised	Single Collimator	0.7 mm Ø
Detector system:	PIN-Diode, 25mm ² with digital pulse processor	Sample stage:	manual, 50 mm height, optional motorised
Radiation safety:	complies with German law for highest safety	Video system	Color CCD-camera with electronic cross hair, field of view 8 x 6 mm, about 20x magnification
Dimensions:	H 35cm, W 40cm, D 60cm	Weight:	about 40 kg
Mains supply:	110V/230V, 50Hz/60Hz		

All specifications are subject to change without notice

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ROENALYTIC GmbH